

Resuscitation Equipment Checklist

- To be reviewed quarterly and after all in office emergencies
- Reviewer to note date of check and initial to confirm completed
- Deficiencies to be recorded on checklist and amended once items replaced

Top Tier Black Trolley

Equipment	Checks	Date	Date	Date	Date
Welch Allyn Vital Sign Machine	Plugged into wall SpO2 cable attached to machine BP cuff cable attached to machine				
BP cuffs	Adult 11 (navy blue) Thigh 13 (black) Large Adult 12 (red) Child 9 (green) Infant (black)				
Pediatric Wireless Pulse Oximeter (green frog)	Turn on to ensure adequate battery life (2xAAA batteries)				
Card with normal ranges for vital signs for reference	None				

Middle Tier Black Trolley

Equipment	Checks	Date	Date	Date	Date
Epinephrine autoinjector* 0.15 mg	Expiry date				
Epinephrine autoinjector* 0.3 mg	Expiry date				
Epinephrine autoinjector* 0.5 mg	Expiry date				
Salbutamol 100mcg x 2	Expiry date				
Rupatadine 10mg tabs	Expiry date Full package (3 tabs)				

Resuscitation Equipment Checklist V1.2 Dr. V. E. Cook

Last update: February 26, 2022

Note: This document is provided as an example. It is recommended that you review the available literature to determine the most appropriate management strategies for your practice.

	At least one box				
Bilastine 20mg tabs	Expiry date Full package (2 tabs) At least one box				
Rupatadine liquid 1mg/mL	Expiry date Full package (30mL) with syringe At least one box				
Bilastine liquid (2.5 mg/mL) – 120 mL bottle	Expiry date 1 box is 120 mL				
Bilastine orodispersible tablets (10 mg)	At least 1 box				
Allerject® trainer	None				
Emerade® trainer	None				
EpiPen® trainer	None				
Aerochamber with mask	In packaging No obvious damage				
Aerochamber with mouthpiece	In packaging No obvious damage				

**Epinephrine autoinjector options include: Allerject® (0.15, 0.3), Emerade® (0.3, 0.5), and EpiPen® (0.15, 0.3)*

Bottom Tier Black Trolley

Equipment	Checks	Date	Date	Date	Date
Pocket mask in white case	Ensure plastic in good quality No damage				
Pediatric non-rebreather mask with O2 tubing (in blue bag)	Unopened No obvious damage				
Pediatric nasal cannula with tubing (in blue bag)	Unopened No obvious damage				
Sure-Flow oxygen tubing (in blue bag)	Unopened No obvious damage				

Sure-Flow oxygen tubing (in black bag)	Unopened No obvious damage				
Adult non-rebreather mask with O2 tubing (in black bag)	Unopened No obvious damage				
Adult oxygen mask with tubing (in black bag)	Unopened No obvious damage				
Set of 8 Guedel Airways in case	40-110mm airways in packaging				
Pediatric Manual Pulmonary Resuscitator (self-inflating bag mask)	Ensure mask, reservoir bag, tubing and 6mm syringe in bag No obvious damage				
Adult Manual Pulmonary Resuscitator (self-inflating bag mask)	Ensure mask, reservoir bag, tubing				
Multi dose inhaler adaptor (blue in Ziploc)	No obvious damage				

Oxygen Tank (located beside fridge in shot room)

Equipment	Check	Date	Date	Date	Date
Oxygen tank	PSI should be between 750-2000 If <750 needs to be refilled				

Prefilled Epinephrine (located in door of fridge in shot room)

Equipment	Check	Date	Date	Date	Date
Prefilled Epinephrine	At least 3x0.5mL with 1" inch needles				

Resuscitation Equipment Checklist V1.2 Dr. V. E. Cook

Last update: February 26, 2022

Note: This document is provided as an example. It is recommended that you review the available literature to determine the most appropriate management strategies for your practice.

syringes with needles	At least 3x<0.5mL(varying increments of 0.05mL) with 5/16"-1.5" needles				
Ampoules of Epinephrine (in top drawer under counter in shot room)	At least 5 ampoules Expiry date Ask doctor to reorder when supply is lower				